

---

## Welcome to the BIOVIA® Interactive Support Matrix

The BIOVIA Interactive Support Matrix is designed to provide a high-level overview of system requirements and does not supersede our published product documentation. This tool can help you answer questions such as:

- When does my BIOVIA product become unsupported? **(01 Support Dates)**
- What client operating systems are compatible with my BIOVIA product? **(02 Client System Requirements)**
- What versions of Microsoft® Office are compatible with my BIOVIA product? **(02 Client System Requirements)**
- What server operating systems are compatible with my BIOVIA product? **(03 Server System Requirements)**
- Do the Application and Database tiers of my BIOVIA product have different operating system requirements? **(03 Server System Requirements)**
- My company is upgrading my desktop Operating System. Will I still be able to run my BIOVIA product? **(04 Client Operating System Support)**
- We are upgrading our Internet browser. Will this affect my use of a BIOVIA product? **(04 Client Operating System Support)**
- When I upgrade my server's operating system, is my BIOVIA product supported on that OS? **(05 Server Operating System Support)**
- Will upgrading to the latest version of Oracle® cause any problems with my BIOVIA product? **(05 Server Operating System Support)**

Your feedback is important to us. If you have any comments or suggestions, please forward it to [BIOVIA Support](#). We want your experience with the Interactive Support Matrix to be everything you need it to be.

## 01 Support Dates

### Browse product release and support dates

This tool displays the release dates, legacy support dates, and unsupported dates for the selected product(s).

There are options to show support dates only for regulated or unregulated environments. You also have the ability to choose to view only the supported versions of a product, all versions, or only the latest version.

The legacy support date is when the product or product version will enter “Legacy Support” as described in the [BIOVIA Software Product Lifecycle and Platform Support Policy](#).

**Support Dates**

Product	Version	Release Date	Legacy Support Beg...	Unsupported as of (...)	Legacy Support Beg...	Unsupported as of (...)
Accelrys Biological Registration	3.1	2/25/2013	12/31/2014	12/31/2015		
Accelrys Chemical Registration	3.0	12/15/2011	5/31/2014	5/31/2015		
Accelrys Chemical Registration	3.1	11/15/2012	5/31/2014	5/31/2015		
Accelrys Chemical Registration	3.2	9/9/2013	9/9/2016	9/9/2017		
Accelrys Draw	3.1 SP1	2/10/2009	2/10/2012	2/10/2013	2/10/2014	2/10/2015
Accelrys Draw	3.2	3/30/2009	3/30/2012	3/30/2013	3/30/2014	3/30/2015
Accelrys Draw	3.2 SP1	6/29/2009	6/29/2012	6/29/2013	6/29/2014	6/29/2015
Accelrys Draw	3.2 SP2	11/13/2009	11/13/2012	11/13/2013	11/13/2014	11/13/2015
Accelrys Draw	3.2 SP3	12/23/2009	12/23/2012	12/23/2013	12/23/2014	12/23/2015
Accelrys Draw	3.3	3/23/2010	3/23/2013	3/23/2014	3/23/2015	3/23/2016
Accelrys Draw	4.0	9/15/2010	9/15/2013	9/15/2014	9/15/2015	9/15/2016
Accelrys Draw	4.0 SP1	4/18/2011	4/18/2014	4/18/2015	4/18/2016	4/18/2017
Accelrys Draw	4.1	7/9/2012	7/9/2015	7/9/2016	7/9/2017	7/9/2018
Accelrys Draw	4.1 SP1	11/9/2012	7/9/2015	7/9/2016	7/9/2017	7/9/2018
Accelrys Electronic Lab Notebook	6.4	12/15/2009	12/31/2012	12/31/2013	12/31/2014	12/31/2015
Accelrys Electronic Lab Notebook	6.4 SP1	2/5/2010	12/31/2012	12/31/2013	12/31/2014	12/31/2015

## 02 Client System Requirements

Browse system requirements for client-based products

This tool will show you client system requirements for BIOVIA® products. This includes compatible operating systems, Internet browsers, and versions of Microsoft Office.

You can select one or more products to view at a time and you may choose to view only the supported versions of a product, all versions, or only the latest version of the selected products.

The screenshot displays the 'Interactive Support Matrix' web application. The left sidebar contains a 'Browse Client System Requirements' section with a product list where 'Accelrys Biological Registration' is selected. Below the list are 'Additional Options' with radio buttons for 'Show All Versions', 'Show Only Latest Version', and 'Show Only Supported Versions'. A 'Submit' button is at the bottom of the sidebar. The main content area, titled 'Client Compatibility', includes a 'Note' at the top, a 'Legend' with color-coded text for legacy and unsupported products, and a table for 'Accelrys Biological Registration'. The table lists version 3.1 with supported operating systems (Windows 7 SP1 and Windows Vista SP2), browsers (Firefox and Internet Explorer), and MS Office requirements (None).

Accelrys Biological Registration			
Version	Client Operating System	Browser(s)	MS Office
3.1	Windows 7 SP1 (32 bit, 64 bit) <i>Windows Vista SP2 (32 bit)</i>	Firefox (32 bit) Internet Explorer 8 (32 bit)	No MS Office Requirements

## 03 Server System Requirements

Browse system requirements for server-based products

This tool will show you server system requirements for BIOVIA® products. This includes compatible operating systems and supported databases.

You can select one or more products to view at a time and you may choose to view only the supported versions of a product, all versions, or only the latest version of the selected products.

The screenshot displays the 'Interactive Support Matrix' web interface. The main content area is titled 'Server Compatibility' and shows requirements for 'Accelrys Chemical Registration'. A legend explains the color coding: orange for legacy support, red for no longer supported, and red with a date for unsupported OS or applications. A table lists supported versions, tiers, operating systems, and databases.

**Note:** See System Requirements documentation for full details.

**Legend:**  
*Accelrys product is in legacy support*  
*Accelrys product is no longer supported*  
*Operating System or third-party application is no longer supported by the vendor*  
 📅 - Support dates are available. Hover mouse over icon for tooltip with the dates.  
 ⓘ - Additional information is available. Hover mouse over icon for tooltip with information.


Accelrys Chemical Registration			
Version	Tier	Server Operating System	Database(s)
3.2 ⓘ	Application	Windows Server 2008 R2 (64 bit) Red Hat Linux Server 5 ES (64 bit)	Oracle 11.2 (64 bit) Oracle 11.2.0 (64 bit)
3.2 ⓘ	Database	Windows Server 2008 R2 (64 bit) Red Hat Linux Server 5 ES (64 bit)	Oracle 11.2 (64 bit) Oracle 11.2.0 (64 bit)
3.1 ⓘ	Application	Windows Server 2008 R2 (64 bit) Red Hat Linux Server 5 ES (64 bit)	Oracle 11.2 (64 bit) Oracle 11.2.0 (64 bit)
3.1 ⓘ	Database	Windows Server 2008 R2 (64 bit) Red Hat Linux Server 5 ES (64 bit)	Oracle 11.2 (64 bit) Oracle 11.2.0 (64 bit)
3.0 ⓘ	Application	Windows Server 2008 R2 (64 bit) Red Hat Linux Server 5 ES (64 bit)	Oracle 11.2 (64 bit) Oracle 11.2.0 (64 bit)
3.0 ⓘ	Database	Windows Server 2008 R2 (64 bit) Red Hat Linux Server 5 ES (64 bit)	Oracle 11.2 (64 bit) Oracle 11.2.0 (64 bit)

**Note:** See System Requirements documentation for full details.

## 04 Client Operating System Support

Browse client operating systems for supported products

This tool lets you choose a client operating system and view all BIOVIA® products that are supported on that operating system.


Interactive Support Matrix
Support Matrix User

**Client Operating System Support**

**Client Operating System:**

Operating System Type  
All

Operating System

- View All Client OS's
- Windows 8.1 (32-bit)
- Windows 8.1 (64-bit)
- SUSE Linux 11 (64-bit)
- Windows 8 (32-bit)
- Windows 8 (64-bit)
- Red Hat Linux 6.0 (64-bit)
- Red Hat Linux 6.3 (64-bit)
- SUSE Linux 10 (64-bit)
- Windows 7 (32-bit)

Submit

**Client OS Compatibility Matrix**

*Note: See System Requirements documentation for full details.*

**Legend:**

- Additional information is available. Hover mouse over icon for tooltip with information.
- Internet Browser Option or Requirement. Hover mouse over for a list of compatible versions.
- MS Office Option or Requirement. Hover mouse over for a list of compatible versions.

Client Operating System Compatibility					
		32-bit			
		Windows			
		8			8.1
		-	Basic	Enterprise	Professional
		-	-	-	-
<a href="#">Accelrys Enterprise Platform</a>	9.0				
	9.1				
<a href="#">Accord for Excel</a>	7.1.5				
<a href="#">BioRails</a>	4.1				
	4.2				
	4.3				
	4.4				
	4.5				
<a href="#">Insight for Excel</a>	2.0				
<a href="#">Materials Studio</a>	7.0				
<a href="#">Morphit</a>	3.0	✓			

*Note: See System Requirements documentation for full details.*

## 05 Server Operating System Support

Browse server operating systems for supported products

This tool lets you choose a server operating system and view all BIOVIA® products that are supported on that operating system.

accelrys® Interactive Support Matrix Beta
Support Matrix User

Home
Parameters
Help

### Server Operating System Support

**Server Operating System:**

Operating System Type  
All

Operating System

- View All Server OS's
- HP UX 11i v3 (64-bit)
- Solaris 10 (32-bit)
- Solaris 10 (64-bit)
- Windows 8 (64-bit)
- Windows 8.1 (64-bit)
- Windows Server 2012 (64-bit)
- SUSE Linux Server 11 (32-bit)
- SUSE Linux Server 11 (64-bit)
- Red Hat Linux Server 6 (64-bit)

Submit

**Note:** See System Requirements documentation for full details.

**Legend:**

- Additional information is available. Hover mouse over icon for tooltip with information.
- Database Option or Requirement. Hover mouse over for a list of compatible database versions.

Server Operating System Compatibility						
		64-bit				
		Windows Server				
		2012				
		-	Datacenter	Essentials	Foundation	Standard
		-	-	-	-	-
<a href="#">Accelrys Enterprise Platform</a>	9.0	-	✓			
	9.1	-	✓			
<a href="#">Discovery Studio</a>	3.5.1	-	✓			
	4.0	-	✓			
<a href="#">Materials Studio</a>	7.0	-		✓	✓	✓

**Note:** See System Requirements documentation for full details.